



Product Code . EL-LG-11168

Thistle Funnel

Description

Thistle Funnel

They are usually passed through a bung or clamped in place.

Thistle Funnel 30 cm stem.




They are particularly useful if the liquid addition must be low down in the flask where a funnel or dropping funnel would have too short a stem.

Thistle Funnels are used to add liquids to a narrow necked flask or other container.

We are leading manufacturers, suppliers of Thistle Funnel for Laboratory Glassware. Contact us to get high quality Thistle Funnel for Laboratory Glassware for schools, colleges, universities, research labs, laboratories and various industries.

```
{ "@context": "https://schema.org/", "@type": "Product", "name": "Thistle Funnel", "image":  
"http://www.educational-equipments.com/images/catalog/product/906132604ThistleFunnel.jpg",  
"description": "They are usually passed through a bung or clamped in place. Thistle Funnel 30 cm  
stem. They are particularly useful if the liquid addition must be low down in the flask where a funnel  
or dropping funnel would have too short a stem. Thistle Funnels are used to add liquids to a narrow  
necked flask or other container. We are leading manufacturers, suppliers of Thistle Funnel for  
Laboratory Glassware. Contact us to get high quality Thistle Funnel for Laboratory Glassware for  
schools, colleges, universities, research labs, laboratories and various industries.", "brand":  
"Educational Lab Equipments", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5",  
"aggregateRating": { "@type": "AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating":
```

"0", "ratingCount": "15" } }

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha, Haryana
Direct Contact Details  +91-98173-19615  sales@educational-equipments.com
 www.educational-equipments.com