

Product Code . EL-BLE-10523

Clinostat Electric



Description

Clinostat Electric

Essential for studies on geotropism and phototropism on seedlings, cut shoots, and small potted plants.

The Clinostat is supplied with a 70mm diameter cork lined/faced disc with a fitting to take a specimen tube.

The angle is adjustable by a thumb screw with a clear plastic cover.

Complete with 3 core cable.

The motor rotates 5 revolutions per hour.




220V, 50Hz. 300mm (H).

The seedlings are pinned to the cork table which is rotated by an electric motor.

Disc diameter: 100mm.

{ "@context": "https://schema.org/", "@type": "Product", "name": "Clinostat Electric", "image": "http://www.educational-equipments.com/images/catalog/product/274435865ClinostatElectric.jpg", "description": "Clinostat Electric Essential for studies on geotropism and phototropism on seedlings, cut shoots, and small potted plants. The Clinostat is supplied with a 70mm diameter cork lined/faced disc with a fitting to take a specimen tube. The angle is adjustable by a thumb screw with a clear plastic cover. Complete with 3 core cable. The motor rotates 5 revolutions per hour. 220V, 50Hz. 300mm (H). The seedlings are pinned to the cork table which is rotated by an electric motor. Disc diameter: 100mm. Clinostat Electric Manufacturers, Clinostat Electric, Clinostat Electric Suppliers,

Clinostat Electric Exporters, Clinostat Electric in India, buy Clinostat Electric, Clinostat Electric Online India, Clinostat Electric Manufacturers in India, Clinostat Electric Biology Lab Equipments", "brand": "Educational Lab Equipments", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating": "0", "ratingCount": "5" }}

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha, Haryana
Direct Contact Details  +91-98173-19615  sales@educational-equipments.com
 www.educational-equipments.com