



Product Code . EL-LG-11069

Adapter Splash Head

Description

Adapter Splash Head

Prevents contents of the flask from being drawn into condenser.

For use with rotary evaporators.

Modified bump traps provide expanded volume fume reaction that have tendency of excessive foam.

Material : Borosilicate Glass.

We are leading manufacturers, suppliers of Adapter Splash Head for Laboratory Glassware. Contact us to get high quality Adapter Splash Head for Laboratory Glassware for schools, colleges, universities, research labs, laboratories and various industries.

```
{ "@context": "https://schema.org/", "@type": "Product", "name": "Adapter Splash Head", "image": "http://www.educational-equipments.com/images/catalog/product/1804936308AdapterSplashHeadWithlogo.jpg", "description": "Adapter Splash Head Prevents contents of the flask from being drawn into condenser. For use with rotary evaporators. Modified bump traps provide expanded volume fume reaction that have tendency of excessive foam. Material : Borosilicate Glass. We are leading manufacturers, suppliers of Adapter Splash Head for Laboratory Glassware. Contact us to get high quality Adapter Splash Head for Laboratory Glassware for schools, colleges, universities, research labs, laboratories and various industries.", "brand": "Educational Lab Equipments", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating": "0", "ratingCount": "10" } }
```

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha, Haryana
Direct Contact Details  +91-98173-19615  sales@educational-equipments.com
 www.educational-equipments.com