



**Product Code . EL-LG-11069**

## Adapter Splash Head

### Description

---

#### Adapter Splash Head

Prevents contents of the flask from being drawn into condenser.

For use with rotary evaporators.

Modified bump traps provide expanded volume fome reaction that have tendency of excessive foam.

Material : Borosilicate Glass.

We are leading manufacturers, suppliers of Adapter Splash Head for Laboratory Glassware. Contact us to get high quality Adapter Splash Head for Laboratory Glassware for schools, colleges, universities, research labs, laboratories and various industries.

```
{ "@context": "https://schema.org/", "@type": "Product", "name": "Adapter Splash Head", "image": "http://www.educational-
```

```
equipments.com/images/catalog/product/1804936308AdapterSplashHeadWithlogo.jpg",  
"description": "Adapter Splash Head Prevents contents of the flask from being drawn into condenser.  
For use with rotary evaporators. Modified bump traps provide expanded volume fome reaction that  
have tendency of excessive foam. Material : Borosilicate Glass. We are leading manufacturers,  
suppliers of Adapter Splash Head for Laboratory Glassware. Contact us to get high quality Adapter  
Splash Head for Laboratory Glassware for schools, colleges, universities, research labs, laboratories  
and various industries.", "brand": "Educational Lab Equipments", "sku": "5", "gtin8": "5", "gtin13": "5",  
"gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating", "ratingValue": "5",  
"bestRating": "5", "worstRating": "0", "ratingCount": "10" } }
```

---

---

Educational Lab Equipments,  
#449, HSIIDC, Industrial Area, Saha, Haryana  
Direct Contact Details  +91-98173-19615  sales@educational-equipments.com  
 [www.educational-equipments.com](http://www.educational-equipments.com)