

**Product Code . EL-EELE-11362**

## Capillary Tube Viscosimeter



### Description

---

#### Capillary Tube Viscosimeter

The liquid whose flow characteristics are to be investigated is placed in a removable thermostatically controlled vessel which may be pressurised. The vessel supplies a capillary tube which is also thermostatically controlled and interchangeable. The apparatus allows Poiseuille flows to be investigated and enables precise measurements to be made of the dynamic viscosities of liquids maintained at constant temperature, using the readings of a single pressure and a recorded time interval. The level of pressurisation of the vessel and the temperature of the liquid can be continuously varied as required ; the liquid is continuously stirred to ensure that its temperature is uniform.

#### Technical Specifications:-




The following instrumentation is provided with the apparatus:

- A stopwatch for flow rate determination with a measuring cylinder.
- Capillary tube diameters : 0.2 and 0.4 mm.
- Temperature range : ambient to 80°C to 0.5°C.
- A thermostatically controlled vessel.
- A mercury manometer.

- A thermometer.
- A magnetic stirrer.
- Pressure range : 0 – 100 cm Hg.
- Dimensions : 37 x 65 x 40 cm high (without manometer) or 120 cm high (with manometer).
- Net weight : 20 kg.

We are leading manufacturers, suppliers of Capillary Tube Viscosimeter for Electronics Engineering Lab Equipments. Contact us to get high quality Capillary Tube Viscosimeter for Electronics Engineering Lab Equipments for schools, colleges, universities, research labs, laboratories and various industries.

```
{ "@context": "https://schema.org/", "@type": "Product", "name": "Capillary Tube Viscosimeter",
"image": "http://www.educational-
equipments.com/images/catalog/product/1380990863CapillaryTubeViscosimeterWithlogo.jpg",
"description": "The liquid whose flow characteristics are to be investigated is placed in a removable
thermostatically controlled vessel which may be pressurised. The vessel supplies a capillary tube
which is also thermostatically controlled and interchangeable. The apparatus allows Poiseuille flows
to be investigated and enables precise measurements to be made of the dynamic viscosities of
liquids maintained at constant temperature, using the readings of a single pressure and a recorded
time interval. The level of pressurisation of the vessel and the temperature of the liquid can be
continuously varied as required ; the liquid is continuously stirred to ensure that its temperature is
uniform. Technical Specifications:- The following instrumentation is provided with the apparatus: • A
stopwatch for flow rate determination with a measuring cylinder. • Capillary tube diameters : 0.2 and
0.4 mm. • Temperature range : ambient to 80°C to 0.5°C. • A thermostatically controlled vessel. • A
mercury manometer. • A thermometer. • A magnetic stirrer. • Pressure range : 0 – 100 cm Hg. •
Dimensions : 37 x 65 x 40 cm high (without manometer) or 120 cm high (with manometer). • Net
weight : 20 kg. We are leading manufacturers, suppliers of Capillary Tube Viscosimeter for
Electronics Engineering Lab Equipments. Contact us to get high quality Capillary Tube Viscosimeter
for Electronics Engineering Lab Equipments for schools, colleges, universities, research labs,
laboratories and various industries.", "brand": "Educational Lab Equipments", "sku": "5", "gtin8": "5",
"gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating",
"ratingValue": "5", "bestRating": "5", "worstRating": "0", "ratingCount": "15" } }
```

Educational Lab Equipments,  
#449, HSIIDC, Industrial Area, Saha, Haryana  
Direct Contact Details  +91-98173-19615  sales@educational-equipments.com  
 www.educational-equipments.com