

**Product Code . EL-TM-11607**

## Pensky Martens Flash Point And Fire Point Apparatus



### Description

---

#### Pensky Martens Flash Point And Fire Point Apparatus




The apparatus consists of Brass Test cup with handle removable cup cover with spring operated rotating shutter having pilot jet, stirrer with flexible shaft. This is widely used for determination of closed cup Flash Point of Fuel Oil, cut back asphalts, other viscous material and suspension of solids having a flash point about 490C (1200F). The assembly rests in Air Bath which is covered with Dome shape metal top. The cup is fitted with insulated handle and locking arrangement near cup plunge. The assembly rests on a round shaped heater with different temperature regulation system suitable for operation on 220 Volts AC mains.

#### Service Required:-

- 15 Amps. Single Phase Electrical Connection.
- Matchbox.
- Ghee.
- Diesel OR 20w/40 Lubricating Oil Sample.
- Suitable Table space to mount apparatus.

We are leading manufacturers, suppliers of Pensky Martens Flash Point And Fire Point Apparatus for Thermodynamics Models. Contact us to get high quality Pensky Martens Flash Point And Fire Point Apparatus for Thermodynamics Models for schools, colleges, universities, research labs, laboratories and various industries.

```
{ "@context": "https://schema.org/", "@type": "Product", "name": "Pensky Martens Flash Point And Fire Point Apparatus", "image": "http://www.educational-equipments.com/images/catalog/product/306913535PenskyMartensFlashPointAndFirePointApparatus.jpg", "description": "The apparatus consists of Brass Test cup with handle removable cup cover with spring operated rotating shutter having pilot jet, stirrer with flexible shaft. This is widely used for determination of closed cup Flash Point of Fuel Oil, cut back asphalts, other viscous material and suspension of solids having a flash point about 490C (1200F). The assembly rests in Air Bath which is covered with Dome shape metal top. The cup is fitted with insulated handle and locking arrangement near cup plunge. The assembly rests on a round shaped heater with different temperature regulation system suitable for operation on 220 Volts AC mains. Service Required:- • 15 Amps. Single Phase Electrical Connection. • Matchbox. • Ghee. • Diesel OR 20w/40 Lubricating Oil Sample. • Suitable Table space to mount apparatus. We are leading manufacturers, suppliers of Pensky Martens Flash Point And Fire Point Apparatus for Thermodynamics Models. Contact us to get high quality Pensky Martens Flash Point And Fire Point Apparatus for Thermodynamics Models for schools, colleges, universities, research labs, laboratories and various industries.", "brand": "Educational Lab Equipments", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating": "0", "ratingCount": "15" } }
```

Educational Lab Equipments,  
#449, HSIIDC, Industrial Area, Saha, Haryana  
Direct Contact Details  +91-98173-19615  sales@educational-equipments.com  
 www.educational-equipments.com