



Product Code . EL-TWL-11760

Chemical Reactor Training Equipment

Description

Chemical Reactor Training Equipment

Description:-

The contact time and reaction temperature can be varied.

Their influence on the conversion of reactants is observed in the respective reactor.

The reaction is moderated using a controlled hot water circuit.

The unit contains two supply tanks for the liquid chemicals.

The properties of different chemical reactors can be investigated by performing a specification reaction using this trainer.

The reactors are available as accessories.

Two variable speed diaphragm pumps pump the reactants through the reactor.




Reactant conversion and reaction temperature are measured with a combined conductivity and temperature measuring unit.

It is also possible to use PC data acquisition.

The resulting reaction product is collected in a stainless steel tray in the unit.

{ "@context": "https://schema.org/", "@type": "Product", "name": "Chemical Reactor Training Equipment", "image": "http://www.educational-equipments.com/images/catalog/product/1352199684

ChemicalReactorTrainingEquipmentWithlogo.jpg", "description": "The contact time and reaction temperature can be varied. Their influence on the conversion of reactants is observed in the respective reactor. The reaction is moderated using a controlled hot water circuit. The unit contains two supply tanks for the liquid chemicals. The properties of different chemical reactors can be investigated by performing a specification reaction using this trainer. The reactors are available as accessories. Two variable speed diaphragm pumps pump the reactants through the reactor. Reactant conversion and reaction temperature are measured with a combined conductivity and temperature measuring unit. It is also possible to use PC data acquisition. The resulting reaction product is collected in a stainless steel tray in the unit.", "brand": "Educational Lab Equipments", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating": "0", "ratingCount": "15" } }

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha, Haryana
Direct Contact Details  +91-98173-19615  sales@educational-equipments.com
 www.educational-equipments.com